

Non sterile

Description	Drug Tariff Code	Unit per Pack
Easi MT 2 litre night bag with 90cm tube, non-return valve, single use twist off drainage outlet	47-60-LBH	10

- Use one per night

Catheter Valves

Valve should be in place for 5-7 days.

Prescribing Advice: 1 box of 5 should be sufficient for 1 month supply.

Unomedical Ltd

Description	Drug Tariff Code	Unit per Pack
Careline Lever Catheter Valve	90-821 1S-20	5
Careline T tap Catheter Valve	90-805 1S-20	5

Penile sheaths

These must be measured and fitted by appropriately trained staff and measuring should only be with manufacturer's measuring guide.

Prescribing Advice: 1 box of 30 should be sufficient for 1 month supply.

HollisterLtd InView Silicone self-adhesive sheath (non latex with anti-kinking funnel)

Description	Drug Tariff Code	Unit per Pack
Special		
25mm	97125	30
29mm	97129	30
32mm	97132	30
36mm	97136	30
41mm	97141	30
Standard		
25mm	97225	30
29mm	97229	30
32mm	97232	30
36mm	97236	30
41mm	97241	30
Extra		
25mm	97325	30
29mm	97329	30
32mm	97332	30
36mm	97336	30
41mm	97341	30

Community Urinary Catheter and Leg Bag and Sheath Formulary



V2

Date of issue March 2019

Date of Review March 2022

Indwelling Catheters

Pre-filled catheters are not recommended for use and any patients currently prescribed these should be switched to one of the foley catheters in the formulary.

Short/medium term use (1-4 weeks) and no latex allergy

Teleflex Medical

RUSCH PTFE AQUAFLATE - PTFE COATED LATEX with sterile water filled syringe for balloon inflation and empty syringe for balloon deflation

Short/medium term use 1-4 weeks (Up to 28 DAYS)

Type and Size	Drug Tariff Code	Unit per Pack
Female 12ch – 10 ml balloon	DP210112	1
Female 14ch – 10ml balloon	DP210114	1
Male 12ch – 10ml balloon	DP310112	1
Male 14ch – 10ml balloon	DP310114	1
Male 16ch – 10ml balloon	DP310116	1

Long term use (3-12weeks)

1st Line Treatment Choice- No Latex Allergy

Teleflex Medical

RUSCH SYMPACATH AQUAFLATE HYDROGEL COATED LATEX with sterile water filled syringe for balloon inflation and empty syringe for balloon deflation. For urethral/supra-pubic use

Type and Size	Drug Tariff Code	Unit per Pack
Female 12ch – 10ml balloon	DH210112	1
Female 14ch – 10ml balloon	DH210114	1
Male 12ch – 10ml balloon	DH310112	1
Male 14ch – 10ml balloon	DH310114	1
Male 16ch – 10ml balloon	DH310116	1

1st Line Treatment Choice- If latex allergy

Teleflex Medical

RUSCH BRILLIANT AQUAFLATE ALL-SILICONE with sterile water filled syringe for balloon inflation and empty syringe for balloon deflation.

Latex free. For urethral/supra-pubic use

Type and Size	Drug Tariff Code	Unit per Pack
All silicone female 12ch – 10 ml balloon	DA210112	1
All silicone female 14ch – 10ml balloon	DA210114	1
All silicone male (standard) 12ch–10ml balloon	DA310112	1
All silicone male (standard) 14ch–10ml balloon	DA310114	1
All silicone male (standard) 16ch–10ml balloon	DA310116	1

Urine Leg Drainage Bags

Bags should be used for 5-7 days

Prescribing advice: 1 box of 10 should be prescribed every 2 months only.

Unomedical Ltd

Careline+ leg bags contain 10 pairs of sterile latex gloves and one pair of soft elasticised cotton Velcro straps per box of 10.

T tap (sterile)

Description	Drug Tariff Code	Unit per Pack
Careline+ 350ml leg bag with short tube t-tap	45-01-SVCG	10
Careline+ 350ml leg bag with long tube t-tap	45-02-LVCG	10
Careline+ 500ml leg bag with short tube t-tap	45-05-SVCG	10
Careline+ 500ml leg bag with long tube t-tap	45-06-LVCG	10
Careline+ 750ml leg bag with short tube t-tap	45-09-SVCG	10
Careline+ 750ml leg bag with long tube t-tap	45-10-LVCG	10

Lever tap (sterile)

Description	Drug Tariff Code	Unit per Pack
Careline+ 500ml leg bag with short tube lever tap	46-05-SVCG	10
Careline+ 500ml leg bag with long tube lever tap	46-06-LVCG	10
Careline+ 750ml leg bag with short tube lever tap	46-09-SVCG	10
Careline+ 750ml leg bag with long tube lever tap	46-10-LVCG	10

Night Drainage bags

Unomedical Ltd

Sterile

Description	Drug Tariff Code	Unit per Pack
Careline+ E4 2 litre night bag with 100cm tube, needle-free sample port, drainable T-tap	45-20-IDCG	10
Careline+ E4 2 litre night bag with 100cm tube, needle-free sample port, drainable lever-tap	46-20-IDCG	10