

### Night Drainage bag: Non sterile – single use bags

Description	Drug Tariff Code	Unit per Pack
Flexicare F2, 2 litre drainage bag with 100cm inlet tube and non-return valve. (single use - tear to open)	00-2202C	10
Flexicare F2 EZ, 2 litre drainage bag with single use snap twist-off tap. (single use - twist-snap tap at top)	00-2203C	10
Flexicare F2 EZB 2 litre drainage bag with single use snap twist-off tap at bottom of bag (single use twist-snap tap at bottom)	00-2208C	10

### Penile sheaths

These must be measured and fitted by appropriately trained staff and measurements should only be with manufacturer's measuring guide.

**Prescribing Advice: 1 box of 30 should be sufficient for 1 month supply.**

Hollister Ltd InView Silicone self-adhesive sheath (non latex with anti-kinking funnel)

Description	Drug Tariff Code	Unit per Pack
<b>Special</b>		
25mm	97125	30
29mm	97129	30
32mm	97132	30
36mm	97136	30
41mm	97141	30
<b>Standard</b>		
25mm	97225	30
29mm	97229	30
32mm	97232	30
36mm	97236	30
41mm	97241	30
<b>Extra</b>		
25mm	97325	30
29mm	97329	30
32mm	97332	30
36mm	97336	30
41mm	97341	30

## Community Urinary Catheter and Leg Bag and Sheath Formulary



**V3**

Date of issue March 2019 (Date amended 18/07/2023)

Date of Review March 2022

## Indwelling Catheters

Pre-filled catheters are not recommended for use and any patients currently prescribed these should be switched to one of the foley catheters in the formulary.

### Short/medium term use (1-4 weeks) and no latex allergy

Teleflex Medical

**RUSCH PTFE AQUAFLATE - PTFE COATED LATEX** with sterile water filled syringe for balloon inflation and empty syringe for balloon deflation

Short/medium term use 1-4 weeks (Up to 28 DAYS)

Type and Size	Drug Tariff Code	Unit per Pack
Female 12ch – 10 ml balloon	DP210112	1
Female 14ch – 10ml balloon	DP210114	1
Male 12ch – 10ml balloon	DP310112	1
Male 14ch – 10ml balloon	DP310114	1
Male 16ch – 10ml balloon	DP310116	1

### Long term use (3-12 weeks)

#### 1<sup>st</sup> Line Treatment Choice- No Latex Allergy

Teleflex Medical

**RUSCH SYMPACATH AQUAFLATE HYDROGEL COATED LATEX** with sterile water filled syringe for balloon inflation and empty syringe for balloon deflation. For urethral/supra-pubic use

Type and Size	Drug Tariff Code	Unit per Pack
Female 12ch – 10ml balloon	DH210112	1
Female 14ch – 10ml balloon	DH210114	1
Male 12ch – 10ml balloon	DH310112	1
Male 14ch – 10ml balloon	DH310114	1
Male 16ch – 10ml balloon	DH310116	1

#### 1<sup>st</sup> Line Treatment Choice- If latex allergy

Teleflex Medical

**RUSCH BRILLIANT AQUAFLATE ALL-SILICONE** with sterile water filled syringe for balloon inflation and empty syringe for balloon deflation.

Latex free. For urethral/supra-pubic use

Type and Size	Drug Tariff Code	Unit per Pack
All silicone female 12ch – 10 ml balloon	DA210112	1
All silicone female 14ch – 10ml balloon	DA210114	1
All silicone male (standard) 12ch–10ml balloon	DA310112	1
All silicone male (standard) 14ch–10ml balloon	DA310114	1
All silicone male (standard) 16ch–10ml balloon	DA310116	1

## Urine Leg Drainage Bags

Bags should be used for 5-7 days

Prescribing advice: 1 box of 10 should be prescribed every 2 months only.

### Lever tap (sterile) (1<sup>st</sup> line option)

Flexicare Medical Ltd

Flexicare leg bag with lever tap outlet, night drainage connector and Velcro anti-slip straps.

Description	Drug Tariff Code	Unit per Pack
Flexicare leg bag with lever tap outlet 350mL, short tube	00-1352	10
Flexicare leg bag with lever tap outlet 500mL, short tube	00-1502	10
Flexicare leg bag with lever tap outlet 750mL, short tube	00-1752	10
Flexicare leg bag with lever tap outlet 350mL, long tube	00-2352	10
Flexicare leg bag with lever tap outlet 500mL, long tube	00-2502	10
Flexicare leg bag with lever tap outlet 750mL, long tube	00-2752	10

### Lever tap or T-tap (sterile) (2<sup>nd</sup> line option)

CliniMed Ltd

CliniSure sterile leg bag with tap, overnight connector and 1 pair of elasticated leg straps.

Description	Drug Tariff Code	Unit per Pack
CliniSure leg bag with 500ml short tube, Lever Tap	CSL500ST	10
CliniSure leg bag with 500ml short tube, T-tap	CST500ST	10
CliniSure leg bag with 500ml long tube, Lever Tap	CSL500LT	10
CliniSure leg bag with 500ml long tube, T-tap	CST500LT	10

### Night Drainage bag: Sterile - can be used for 5-7 days

Description	Drug Tariff Code	Unit per Pack
Flexicare F4, Sterile 2 litre drainage bag with 100cm inlet tube, non-return valve and tap outlet.	00-1200	10
Flexicare F4L, Sterile 2 litre drainage bag with 100cm inlet tube, non-return valve and lever tap outlet.	00-1201	10